



# Bias Cable from detector fanout board to cryostat connector



Micro D-25 way Cable mounted socket

16-26P Cryostat mounted plug



## Notes

- Cable length = 35 cm
- Manganin Wiring
- Use Teflon heat shrink if required

Title		
Internal ESO cryostat wiring for CCD201 bias/video cable		
Size	Number	Revision
A3		Version 1
Date:	21/01/2008	Sheet of
File:	U:\L3CCD\ESO_CRYOSTAT_WIRING	DRASEBY.V1.SCHDOC