


# Cable from CCD201 Detector fanout board to socket on cryostat wall

Micro D-37 way Cable mounted socket

20-41P Amphenol cryostat mounted plug



**Notes**

**USE CONSTANTAN WIRE**

Cable is internal to Cryostat  
Cable length = 35 cm  
Use Teflon heat shrink only

Modified to mark up twisted pairs

Title		
Internal ESO cryostat wiring for CCD201 clock cable		
Size	Number	Revision
A3		Version 1
Date:	21/01/2008	Sheet of
File:	U:\L3CCD\...ESO_Cryostat_Wiring_Clocks\Drawings	1 of 1